

Nome do usuário:

Data

Peso:

Altura:

Diagnóstico:

Informações Adicionais

Option #	HCPCS	Definition	Price	Option #	HCPCS	Definition	Price
----------	-------	------------	-------	----------	-------	------------	-------

1. CONFIGURAÇÃO ESTRUTURA

Modelo base

<input checked="" type="checkbox"/>	186WM60	K0005	Quickie Q7 - Adjustable (265 lbs)	\$
-------------------------------------	---------	-------	-----------------------------------	----

Kit Performance

<input checked="" type="checkbox"/>	186QD01	Garfo Unilateral	\$
		Roda ultraleve	\$
		Pneu Schwalbe Marathon Plus	\$
		Freio Compacto	\$

Kit Ultralite

<input checked="" type="checkbox"/>	186QD02	Lona assento Ultralite	\$
		Roda ultraleve	\$
		Pneu Schwalbe Speed Run Tire	\$
		Freio Compacto	\$

Angulação Estrutura

See **FIGURE 1** Measure from front frame to floor.

<input type="checkbox"/>	186FA1	75°	NC
<input checked="" type="checkbox"/>	186FA8	80°	NC
<input type="checkbox"/>	186FA2	85°	NC

Transito

<input checked="" type="checkbox"/>	186TO45	Kit Transito	\$
-------------------------------------	---------	--------------	----

Fechamento Frontal

See **FIGURE 2** Frame Inset is the distance from the outside of the seat tube to the outside of the footrest.

<input type="checkbox"/>	186FF1	0"	NC
	186FF2	1" (NA w/ 12" wide)	NC
<input type="checkbox"/>	186FF3	2" (NA with 12" & 13" wide)	NC

Estilo Estrutura

See **FIGURE 3** Measurement for set up.

<input type="checkbox"/>	186FS9	Standard	NC
<input type="checkbox"/>	186FS34	Classico	NC

* **Standard** – A estrutura é 4cm menor que a clássica **Clássico – Tradicional** posicionamento utilizado para a distância entre o assento e a curva frontal. A profundidade do assento não é alterada para ambos os modelos.

Posição da Roda Dianteira

See **FIGURE 4**

<input type="checkbox"/>	186LH01	Setting 1	Sugerida para rodas 3" & 4"	NC
<input checked="" type="checkbox"/>	186LH02	Setting 2	Sugerida para rodas 5" & 6"	NC
	186LH03	Setting 3	NA w/ 6" casters	NC

Direct Mount (Opcional)

Apenas válido com freio compacto e cambagem de 0°.

<input type="checkbox"/>	186WH13	Direct Mount Wheel Lock	NC
--------------------------	---------	-------------------------	----

Figure 1




Figure 2




Figure 3





Figure 4



1. CONFIGURAÇÃO ESTRUTURA

Encosto ajustável em profundidade (Opcional)

See **Figure 5** Depth Adjust Backrest extends the frame rear of the seat sling, select a factory setting.

<input type="checkbox"/>	186DP02	0 Setting Backrest at standard position, field adjustable to -1" or -2"	\$165.00
<input type="checkbox"/>	186DP03	-1 Setting Backrest -1" offset, field adjustable to flush or -2"	\$165.00
<input type="checkbox"/>	186DP04	-2 Setting Backrest -2" offset, field adjustable to flush or -1"	\$165.00




Figure 5

0 Position

-1" Position

-2" Position

Largura do assento

See **FIGURE 6** Measure from outside of seat tube to the outside of opposite seat tube.


<input type="checkbox"/>	186S3	30cm	NC	<input type="checkbox"/>	186S8	44cm	NC
<input type="checkbox"/>	186S4	33cm	NC	<input type="checkbox"/>	186S9	46cm	NC
<input type="checkbox"/>	186S5	36cm	NC	<input type="checkbox"/>	186S10	48cm	NC
<input type="checkbox"/>	186S6	38cm	NC	<input type="checkbox"/>	186S48	51cm (E2201)	\$
<input checked="" type="checkbox"/>	186S7	40cm	NC				

Profundidade do assento

See **FIGURE 7** Measure from front of back post to front edge of seat sling. Seat sling is 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.

<input type="checkbox"/>	186SD5	30cm	NC	<input type="checkbox"/>	186SD10	44cm	NC
<input checked="" type="checkbox"/>	186SD6	33cm	NC	<input checked="" type="checkbox"/>	186SD11	46cm	NC
<input checked="" type="checkbox"/>	186SD7	36cm	NC	<input checked="" type="checkbox"/>	186SD12	48cm	NC
<input checked="" type="checkbox"/>	186SD8	38cm	NC	<input checked="" type="checkbox"/>	186SD13	51cm (E2203)	\$
<input checked="" type="checkbox"/>	186SD9	40cm	NC				

Figure 6




Frame Inset

Seat Width

Footrest Width

Figure 7



Seat Depth

Frame Angle

Opções assento

<input checked="" type="checkbox"/>	186OP01	Lona Standard	NC
<input checked="" type="checkbox"/>	186OP02	Lona Ultra fina Válido de 36-44cm profundidade e larguar.	\$
<input checked="" type="checkbox"/>	186OP04	Omit Seat Sling	NC

Amofada Jay (Opcional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

<input checked="" type="checkbox"/>	186JC20	JAY Bundled Cushion	\$
		Modelo:	
		Largura:	Profundidade:

2. Altura Frontal

Altura frontal

See **FIGURE 8** Measure from the front of the seat rail to the floor.

<input type="checkbox"/>	186SH14	40cm	NC	<input type="checkbox"/>	186SH20	47cm	NC
<input checked="" type="checkbox"/>	186SH15	42cm	NC	<input checked="" type="checkbox"/>	186SH21	48cm	NC
<input checked="" type="checkbox"/>	186SH16	43cm	NC	<input checked="" type="checkbox"/>	186SH23	49cm	NC
<input checked="" type="checkbox"/>	186SH17	44cm	NC	<input checked="" type="checkbox"/>	186SH24	51cm	NC
<input checked="" type="checkbox"/>	186SH18	46cm	NC	<input checked="" type="checkbox"/>	186SH26	52cm	NC
				<input checked="" type="checkbox"/>	186SH27	53cm	NC

Figure 8



Rear Seat Height

Front Seat Height

Rodas Dianteiras

<input checked="" type="checkbox"/>	186CW1	3" Micro	NC
<input checked="" type="checkbox"/>	186CW3	4" Low Profile Polyurethane	NC
<input checked="" type="checkbox"/>	186CW31	4" Micro Lighted	\$
<input checked="" type="checkbox"/>	186CW72	4" x 1" Aluminum Poly	\$
<input checked="" type="checkbox"/>	186CW44	4" x 1.25 Semi Pneumatic	\$
<input checked="" type="checkbox"/>	186CW64	4" x 1.5" Aluminum Soft Roll	\$
<input checked="" type="checkbox"/>	186CW5	5" Low Profile Polyurethane	NC
<input checked="" type="checkbox"/>	186CW51	5" Micro Lighted	\$
<input checked="" type="checkbox"/>	186CW6	5" Soft Roll (Grey) -NA w/Frog Legs	\$
<input checked="" type="checkbox"/>	186CW65	5" x 1" Aluminum Poly	\$
<input checked="" type="checkbox"/>	186CW7	5" x 1.5" Semi Pneumatic	\$
<input checked="" type="checkbox"/>	186CW66	5" x 1.5" Aluminum Soft Roll	\$
<input checked="" type="checkbox"/>	186CW34	6" Poly	NC
<input checked="" type="checkbox"/>	186CW33	E2214 E2215 6" Pneumatic	\$
<input checked="" type="checkbox"/>	186CW35	6" x 1.5" Semi Pneumatic	\$
<input checked="" type="checkbox"/>	186CW68	6" x 1.5" Aluminum Soft Roll	\$

Garfos e Suspensões

<input checked="" type="checkbox"/>	186CF74	Carbotecture Forks	NC
<input checked="" type="checkbox"/>	186CF65	Múltiplas posições NEW!	NC
<input checked="" type="checkbox"/>	186CF69	Múltiplas posições - unilateral	\$
<input checked="" type="checkbox"/>	186CF53	E1015 Frog Legs 0 até 45kg.	\$
<input checked="" type="checkbox"/>	186CF54	E1015 Frog Legs 46kg- 86kg	\$
<input checked="" type="checkbox"/>	186CF55	E1015 Frog Legs 87kg-120kg.	\$

3. Apoio de pés

Apoio dos pés				
<input checked="" type="checkbox"/>	186F229		Alumínio Tubular	NC
<input checked="" type="checkbox"/>	186F274		Alumínio Tubular com base	\$
<input checked="" type="checkbox"/>	186F276	K0040	Ângulo Ajustável	\$
<input checked="" type="checkbox"/>	186F51	K0040	Pedal Plataforma	\$
<input checked="" type="checkbox"/>	186F111	K0040	Pedal Plataforma Ultralite	\$

Opcionais				
<input checked="" type="checkbox"/>	186DI08	K0037	High Mount Ext Tubes (NA w/ Tubular footrest)	\$
<input checked="" type="checkbox"/>	186DI03	K0038	Faixa para perna	\$
<input checked="" type="checkbox"/>	186DI04	K0038	Faixa extra para perna	\$
<input checked="" type="checkbox"/>	186DI05		Capa em gel para hanger	\$
<input checked="" type="checkbox"/>	186DI07		Neoprene Impact Guards	\$
<input checked="" type="checkbox"/>	186DI02		Wear Pads	\$
<input checked="" type="checkbox"/>	186DI06		Position Plates (NA for Tubular Footrest)	\$

NEW!

Adaptive Platform Spacer

NA with Aluminum Tubular or Ultra Lite Platform. Now available with Angle Adj. FlipBack.

<input checked="" type="checkbox"/>	186AD1		Adaptive Platform 2" Spacer	\$
<input checked="" type="checkbox"/>	186AD2		Adaptive Platform 4" Spacer	\$
<input type="checkbox"/>	186AD3		Adaptive Platform 6" Spacer	\$

4. Encosto

Back Type				
<input checked="" type="checkbox"/>	186BT96		Dobrável ângulo ajustável	NC
<input checked="" type="checkbox"/>	186BT3		Não dobrável ângulo Ajustável	NC

Choose Option 1 or 2

1. Encosto Altura Ajustável				2. Encosto Fixo Ultralight									
<input checked="" type="checkbox"/>	186B110	Setting	25-36cm Altura	NC	Select back height. 8", 9", or 10" only available with omit push handles. 17", 18", or 19" only available with fold down push handles								
<input checked="" type="checkbox"/>	186B110	Setting	30-40cm Altura	NC				<input type="checkbox"/>	186BB21	8"	<input type="checkbox"/>	186BB17	36
<input checked="" type="checkbox"/>	186B110	Setting	36-46cm Altura	NC				<input type="checkbox"/>	186BB20	9"	<input type="checkbox"/>	186BB18	38
<input checked="" type="checkbox"/>	186B110	Setting	40-51cm Altura	NC				<input type="checkbox"/>	186BB13	10"	<input type="checkbox"/>	186BB19	40
<input checked="" type="checkbox"/>	186B110	Setting		NC				<input type="checkbox"/>	186BB14	11"	<input type="checkbox"/>	186BB10	44
					<input type="checkbox"/>	186BB15	12"	<input type="checkbox"/>	186BB11	46			
					<input type="checkbox"/>	186BB16	13"	<input type="checkbox"/>	186BB12	19"			

Barra				
<input checked="" type="checkbox"/>	186QT01		Standard	NC
<input checked="" type="checkbox"/>	186QT02		Low Profile	\$

Manoplas				
<input checked="" type="checkbox"/>	186PH6		Integral	NC
<input checked="" type="checkbox"/>	186PH5		Retrátil (dobrável)	\$
<input type="checkbox"/>	186PH2		Omitir manoplas	NC

Lona Encosto				
<input checked="" type="checkbox"/>	186BU1		Standard	NC
<input checked="" type="checkbox"/>	186BU70	E2611	Tensão Ajustável Nylon Balístico	\$
<input checked="" type="checkbox"/>	186BU71	E2611	Tensão Ajustável 3DX Vented	\$
<input type="checkbox"/>	186BU14		Omitir Lona	NC

Encosto J3				
<input checked="" type="checkbox"/>	186JB01		J3 Back	\$

Broches militares				
Only available on tension adjustable upholstery and Jay 3 Back				
<input type="checkbox"/>	186VU1		US Marines Patch	\$
<input type="checkbox"/>	186VU2		US Army Patch	\$
<input type="checkbox"/>	186VU3		US Navy Patch	\$
<input type="checkbox"/>	186VU4		US Air Force Patch	\$
<input type="checkbox"/>	186VU12		Coast Guard Patch	\$
<input type="checkbox"/>	186VU5		Purple Heart Patch	\$
<input type="checkbox"/>	186VU8		Purple Heart with Marines Patch	\$
<input type="checkbox"/>	186VU9		Purple Heart with Army Patch	\$
<input type="checkbox"/>	186VU10		Purple Heart with Navy Patch	\$
<input type="checkbox"/>	186VU11		Purple Heart with Air Force Patch	\$
<input type="checkbox"/>	186VU13		Purple Heart with Coast Guard Patch	\$

NEW!

Opcionais para Encosto

Telescoping Push Handles have Double "D" tubing to eliminate rotation.

<input type="checkbox"/>	186OD5		Narrow Backrest Inset	\$
<input checked="" type="checkbox"/>	186A192		Manopla Telescópica - Baixa	\$
<input checked="" type="checkbox"/>	186A195		Manopla Telescópica - Alta	\$

5. Configuração Traseira

Eixos Traseiros

<input type="checkbox"/>		Centro de Gravidade Ajustável (1"-4")	NC
<input type="checkbox"/>		Centro de Gravidade Ajustável (0"-3")	NC

Suspensão Traseira (Opcional)

<input type="checkbox"/>				
<input type="checkbox"/>	186SV08	Weight from 0 lbs to 175 lbs		\$
<input type="checkbox"/>	186SV09	Weight from 176 lbs. to 265 lbs		\$

Xtender Power Assist

Rear seat height is 16.5" to 20", 0° camber, ships with 1" COG, pneumatic tire. 24" wheels only. Anti-Tips no charge. NA w/ Padded Swing Arms or Transit.

NEW!

<input type="checkbox"/>	186XT2	E0986	Xtender Power Assist 4X Multiplier NICKEL	\$
<input checked="" type="checkbox"/>	186XT3	E0986	Xtender Power Assist 4X Multiplier LITHIUM	\$

Please select one of the following handrim options

<input type="checkbox"/>	186XH3		Plastic Coated	NC
<input type="checkbox"/>	186XH4		Stainless Steel	NC

Xtender – Outras Opções

NEW!

<input checked="" type="checkbox"/>	186XA1		Spare NiMH Power Pack - NICKEL	\$
<input checked="" type="checkbox"/>	186XA7		Spare NiMH Power Pack - LITHIUM	\$

Xtender Opções de Rodas

NEW!

<input checked="" type="checkbox"/>	186XA2	E2211 E2212	SPARE 24" Lite Spoke Pneu Alum	NC
<input checked="" type="checkbox"/>	186XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic	\$

Posicionamento do Centro de gravidade

Measure from front of back post to the center of rear axle.

<input type="checkbox"/>	186CC12	0		NC
<input type="checkbox"/>	186CC1	1cm		NC
<input type="checkbox"/>	186CC2	2,5cm		NC
<input type="checkbox"/>	186CC3	4cm		NC
<input type="checkbox"/>	186CC4	5 cm		NC
<input type="checkbox"/>	186CC5	6cm		NC
<input type="checkbox"/>	186CC6	8cm		NC
<input type="checkbox"/>	186CC7	9cm		NC
<input type="checkbox"/>	186CC8	10cm		NC

Cambagem

<input checked="" type="checkbox"/>	186CI1		0°	NC
<input checked="" type="checkbox"/>	186CI2		3°	NC
<input checked="" type="checkbox"/>	186CI4		6°	NC

Wheel	5 Spoke Mag	Mag	Lite Spoke	Ultra Lt wt	Spinergy
Full Poly Tire	20,22,24	20,22,24	20,22,24		
High Press. Clincher			22,24,25,26	22,24,25	22,24,25,26
SHOX Perform. Poly	20,22,24	24	24,25		24,25
Knobby V-TRAK			24		24
Pneumatic	20,22,24	20,22,24	20,22,24,26		22,24,26
Pneumatic w/Airless Insert	20,22,24	20,22,24	20,22,24,26		22,24,26
Schwalbe Marathon			24,25,26	22,24,25	24,25,26
Schwalbe Speed Run			24,25	24,25	24,25
22" = 489 mm 24" = 540 mm 25" = 559 26" = 590					

Matrix altura Anterior x Posterior

30cm Profundidade	2.5" Max Front to Rear Seat Hgt.
33,36cm Profundidade	3.0" Max Front to Rear Seat Hgt.
38,40 cm Profundidade	3.5" Max Front to Rear Seat Hgt.
44,46 cm Profundidade	4.0" Max Front to Rear Seat Hgt.
48,51 cm Profundidade	4.5" Max Front to Rear Seat Hgt.

Altura Posterior

<input checked="" type="checkbox"/>	186RH11		33 – Somente com rodas traseiras de 20"	NC
<input checked="" type="checkbox"/>	186RH12		34 – Somente com rodas traseiras de 20"	NC
<input checked="" type="checkbox"/>	186RH13		36 _ Com rodas tras. De 20"	NC
<input checked="" type="checkbox"/>	186RH14		37 - . Com rodas tras. De 20"	NC
<input checked="" type="checkbox"/>	186RH15		38 – Somente com rodas traseirasde 20 e 22"	NC
<input checked="" type="checkbox"/>	186RH16		39 - Somente com rodas traseirasde 20 e 22"	NC
<input checked="" type="checkbox"/>	186RH17	40	ND com rodas 24, 25 & 26"	NC
<input checked="" type="checkbox"/>	186RH18	41	ND com rodas 24, 25 & 26"	NC
<input checked="" type="checkbox"/>	186RH19	44		NC
<input checked="" type="checkbox"/>	186RH26	45		NC
<input checked="" type="checkbox"/>	186RH20	46		NC
<input checked="" type="checkbox"/>	186RH27	47		NC
<input checked="" type="checkbox"/>	186RH21	48		NC
<input checked="" type="checkbox"/>	186RH37	49		NC
<input checked="" type="checkbox"/>	186RH22	51	ND com suspensão	NC

Axle

<input checked="" type="checkbox"/>	186AX1		Eixo de liberação rápida	NC
<input checked="" type="checkbox"/>	186AX13		Eixo de liberação rápida nas 4 rodas	\$

Rodas Traseiras

Ultra Lightweight apenas disponível com aro de propulsão em alumínio

<input type="checkbox"/>	186RW3		20" Mag	NEW!	NC
<input checked="" type="checkbox"/>	186RW8		22" Mag	NEW!	NC
<input checked="" type="checkbox"/>	186RW1		24" Mag	NEW!	NC
<input checked="" type="checkbox"/>	186RW103		20" Mag Wheel 5 Spoke		NC
<input checked="" type="checkbox"/>	186RW104		22" Mag Wheel 5 Spoke		NC
<input checked="" type="checkbox"/>	186RW102		24" Mag Wheel 5 Spoke		NC
<input checked="" type="checkbox"/>	186RW90		20" Lite Spoke		\$
<input checked="" type="checkbox"/>	186RW91		22" Lite Spoke		\$
<input checked="" type="checkbox"/>	186RW92		24" Lite Spoke		NC
<input checked="" type="checkbox"/>	186RW93		25" Lite Spoke		NC
<input checked="" type="checkbox"/>	186RW94		26" Lite Spoke		\$
<input checked="" type="checkbox"/>	186RW06		22" Ultra Lightweight Spoke 2.0		\$
<input checked="" type="checkbox"/>	186RW10		24" Ultra Lightweight Spoke 2.0		\$
<input checked="" type="checkbox"/>	186RW16		25" Ultra Lightweight Spoke 2.0		\$
<input checked="" type="checkbox"/>	186RW28		22" Spinergy SPOX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW95		24" Spinergy SPOX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW96		25" Spinergy SPOX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW97		26" Spinergy SPOX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW32		22" Spinergy LX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW98		24" Spinergy LX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW99		25" Spinergy LX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW40		26" Spinergy LX Black Spokes		\$
<input checked="" type="checkbox"/>	186RW50		24" Mountain		\$

5. CONFIGURAÇÃO TRASEIRA

Pneu Traseiro				
<input checked="" type="checkbox"/>	186RT3	E2218	Full Poly Tire	NC
<input checked="" type="checkbox"/>	186RT6	E2211 E2212	High-Pressure Clincher	\$
<input checked="" type="checkbox"/>	186RT38	E2218	SHOX Performance Poly	\$
<input checked="" type="checkbox"/>	186RT27	E2211 E2212	Knobby V-TRAK by Kenda	\$
<input checked="" type="checkbox"/>	186RT12	E2211 E2212	Mountain <i>Only available in 24"</i>	\$
<input checked="" type="checkbox"/>	186RT1	E2211 E2212	Pneumatic	NC
<input checked="" type="checkbox"/>	186RT2	E2211 E2213	Pneumatic Airless Insert	\$
<input checked="" type="checkbox"/>	186RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$
<input checked="" type="checkbox"/>	186RT28	E2211 E2212	Schwalbe Speed Run Tire	\$

Aro de Propulsão				
<i>Natural Fit Handrims compatível com rodas 22", 24", & 25".</i>				
<input type="checkbox"/>	186HR1		Aluminum Anodized	NC
<input type="checkbox"/>	186HR4		Plastic Coated ND c/ rodas 25" & 26" Lite Spoke	\$
<input checked="" type="checkbox"/>	186HR22		Natural Fit - Standard Grip	\$
<input checked="" type="checkbox"/>	186HR23		Natural Fit - Super Grip	\$
<input checked="" type="checkbox"/>	186HR30		Natural Fit LT - Standard Grip	\$
<input checked="" type="checkbox"/>	186HR31		Natural Fit LT - Super Grip	\$
<input checked="" type="checkbox"/>	186HR24		Spinergy Flex Rim <i>Apenas disponível c/ Spinergy LX 24" & 25"</i>	\$
<input checked="" type="checkbox"/>	186HR7		Omit Handrim <i>ND c/ rodas Ultra Lightweight</i>	NC

5. APOIO DE BRAÇOS

Apoio de braço acolchoado- Swing Away							
<i>Some limitations when used with the Angle Adjustable Back in the rear most setting</i>							
<i>Please select a setting, will default to 11" if no selection is made</i>							
<input type="checkbox"/>	186AA18	E0973	Acolchoado Swing-Away STD	\$			
	Setting	<input type="checkbox"/> 23	<input type="checkbox"/> 25	<input type="checkbox"/> 28	<input type="checkbox"/> 30	<input type="checkbox"/> 33	<input type="checkbox"/> 36

Armpad para Apoio altura ajustável				
<input checked="" type="checkbox"/>	186AM23		Desk Length (25cm) <i>Classic</i>	NC
<input checked="" type="checkbox"/>	186AM21		Desk Length (25cm) <i>Contour</i>	NC
<input checked="" type="checkbox"/>	186AM24		Full Length (36cm) <i>Classic</i>	NC
<input checked="" type="checkbox"/>	186AM22		Full Length (36cm) <i>Contour</i>	NC

Apoio de braço - altura ajustável				
<input checked="" type="checkbox"/>	186AA11	E0973	Single Post Height Adjustable - Ajuste de 20 -33cm altura	\$
<input checked="" type="checkbox"/>	186AA29	E0973	Single Post Height Adjustable - Low Ajuste de 16- 26,5cm altura	\$

Opcionais				
<input checked="" type="checkbox"/>	186NP01		Pad de gel para apoio swing - away	\$

6. ACESSÓRIOS

Freios				
<i>Not all wheel locks available in all seat depths, wheel size, or COG</i>				
<input type="checkbox"/>	186WL1		Freio alto acion. Para frente	NC
<input checked="" type="checkbox"/>	186WL2		Freio alto acion. Para trás	NC
<input checked="" type="checkbox"/>	186WL45		Ergo Scissor	NC
<input checked="" type="checkbox"/>	186WL30		Compacto	\$
<input checked="" type="checkbox"/>	186WL10		Omit Wheel Lock	NC

Protetor de roupas - Cor				
<i>Available with Aluminum or Aluminum with Fender.</i>				
<input type="checkbox"/>	186GE01		Preto	NC
<input checked="" type="checkbox"/>	186GE02		Mesma cor da cadeira - cores padrão	NC
<input checked="" type="checkbox"/>	186GE03		Mesma cor da cadeira - cor customiz.	\$

Opcionais				
<input checked="" type="checkbox"/>	186WL12		Sem montagem	NC

Cintos de posicionamento				
<input checked="" type="checkbox"/>	186A117	E0978	Belt 2" Aircraft Buckle w/clamps	\$
<input checked="" type="checkbox"/>	186A118	E0978	Belt 2" Aircraft Buckle Padded w/clamps	\$
<input checked="" type="checkbox"/>	186A141	E0978	Belt 1.5" Auto Buckle Universal w/clamps	\$

Extensor de freio				
<i>Ext. handles only available with high mount push or pull wheel locks</i>				
<input checked="" type="checkbox"/>	186EH7	E0961	Extension Handle (Pair)	\$65.00

Protetor de Rodas				
<i>Not available with 5 spoke MAG wheels</i>				
	186A23	K0065	Preto	\$
<input checked="" type="checkbox"/>	186A302	K0065	Quickie "Q" Clear	\$

Anti-tombos				
<input checked="" type="checkbox"/>	186AT37	E0971	Anti Tips (Pair)	\$

Acessórios				
<input checked="" type="checkbox"/>	186A1		Mochila adulto	\$
<input checked="" type="checkbox"/>	186A2		Mochila infantil	\$
<input checked="" type="checkbox"/>	186A8		Black Seat Pouch	\$45.00
<input checked="" type="checkbox"/>	186A203	E2207	Suporte muletas	\$
<input checked="" type="checkbox"/>	186A83		Kit ferramentas	\$

Protetor de Roupas				
<input type="checkbox"/>	186SG33	K0108	Plástico Standard <i>ND com apoio de braço Ang Adjs</i>	\$
<input checked="" type="checkbox"/>	186SG29	K0108	Plástico Low <i>ND com apoio de braço</i>	\$
<input checked="" type="checkbox"/>	186SG35	K0108	Alumínio <i>ND com apoio de braço Ang Adjs, freio alto</i>	\$
<input checked="" type="checkbox"/>	186SG36	K0108	Alumínio com fenda, <i>ND com apoio de braço Ana Adjs e freio alto</i>	\$
<input checked="" type="checkbox"/>	186SG12	K0108	Fibra de carbon, <i>ND com apoio de braço Ang Adjs e freio alto</i>	\$



6. ACESSÓRIOS

Cor da Estrutura			
<input type="checkbox"/>	186C118	Bead Blast	NC
<input type="checkbox"/>	186C140	Aztec Gold	NC
<input type="checkbox"/>	186C2	Black <i>(Not avail as base w/Zebra)</i>	NC
<input type="checkbox"/>	186C141	Black Cherry	NC
<input type="checkbox"/>	186C10	Black Opal	NC
<input type="checkbox"/>	186C38	Blue Opal	NC
<input type="checkbox"/>	186C36	Candy Blue	NC
<input type="checkbox"/>	186C14	Candy Purple	NC
<input type="checkbox"/>	186C35	Candy Red	NC
<input type="checkbox"/>	186C142	Electric Blue	NC
<input type="checkbox"/>	186C143	Evergreen	NC
<input type="checkbox"/>	186C126	Glow	NC
<input type="checkbox"/>	186C144	Green Apple	NC
<input type="checkbox"/>	186C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	186C148	Mauve Pink	NC
<input type="checkbox"/>	186C103	Rootbeer	NC
<input type="checkbox"/>	186C114	Sparkle Silver	NC
<input type="checkbox"/>	186C145	Sunrise Orange	NC
<input type="checkbox"/>	186C113	Titanium Color Paint	NC
<input type="checkbox"/>	186C8	Yellow	NC
<input type="checkbox"/>	186C44	Glossy White - Limited Edition	NC

Cores Foscas			
<input type="checkbox"/>	186C116	Matte Black	NC
<input type="checkbox"/>	186C151	Matte Black Cherry	NC
<input type="checkbox"/>	186C150	Matte Electric Blue	NC
<input type="checkbox"/>	186C149	Matte Evergreen	NC
<input type="checkbox"/>	186C152	Matte Purple	NC

Cores customizadas			
<input type="checkbox"/>	186C112	Desert Camouflage	\$
<input type="checkbox"/>	186C43	Mossy Oak Camouflage	\$
<input type="checkbox"/>	186C42	Stars & Stripes	\$
<input type="checkbox"/>	186C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Black</i>	\$

Peças anodizadas			
Garfos, Suporte do Encosto e abraçadeira do suporte de roda			
<input type="checkbox"/>	186EB07	Preto	NC
<input type="checkbox"/>	186EB08	Azul	\$
<input type="checkbox"/>	186EB09	Verde	\$
<input type="checkbox"/>	186EB10	Roxo	\$
<input type="checkbox"/>	186EB11	Vermelho	\$
<input type="checkbox"/>	186EB12	Cinza Titânio	\$

Built For Me

*To inquire about any unique customization for Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 3339011*



Commonly Requested Options Include:

- ✓ Swing Away Knee Adductors
- ✓ One Arm Drive
- ✓ Unilateral Wheel Locks
- ✓ Seat Pan
- ✓ Caster Pin Lock
- ✓ Grade Aids
- ✓ Spinergy Color Spokes
- ✓ Double D Telescoping Push Handles

